Serafol[®] ABO

Serafol[®] ABO+D

Bedside-cards for identity proof (ready to use twin-cards)

Information for use by trained personnel

Test purpose

Serafol[®] ABO and Serafol[®] ABO+D are used for an identity proof immediately prior to a blood transfusion (bedside test). The test is the confirmation of earlier ABO- and D-blood typing of the recipient and ensures the compatibility of blood types of the recipient and the blood to be transfused. Thus possible mismatches can be detected.

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Test principle

Test cards coated with dried specific monoclonal reagents.

The test principle is a hemagglutination test for the detection of the respective red cell antigens. ABO blood types and the rhesus characteristic D are defined by the presence or absence of the antigens A, B, and D on red cells. If such antigens are present on the red cells, they will be agglutinated by the corresponding antibodies (positive reaction).

Twin card / Single card

Serafol cards can be used either for ABO, ABO+D identity tests of recipients and blood donation (twin card) or for tests of 2 recipients (single card).

Composition

Serafol[®] ABO: The single card contains 3 reaction fields: one coated with Anti-A (clone A003, ≥ 200 IU/ml) and one coated with Anti-B (clone B005, ≥ 200 IU/ml) as well as one field for the blood to be tested or auto control.

Serafol[®] ABO+D: The single card contains 4 reaction fields: one coated with Anti-A (clone A003, ≥ 200 IU/ml), one coated with Anti-B (clone B005, ≥ 200 IU/ml), and one coated with Anti-D (clone BS226, Titre: $\geq 1:32$) as well as one field for the blood to be tested or auto control.

Stabilizer: ≤ 0.1 % NaN₃

The antibodies are derived from culture supernatants of stable cell lines and demonstrate the consistent specificity and reproducibility characteristic for monoclonal antibodies.

Shelf life

Under proper storage conditions Serafol cards wrapped in foil can be used until the expiry date printed both on the foil and the outer label.

Serafol cards should be taken from their foil-wrapper only immediately before use. Unused cards that were taken from the foil or unused halves of twin cards should be placed back in the foil and be used within 7 days.

Additionally needed reagents and materials

Isotonic saline solution or water

Pipettes

	Catalogue no.	Package size
Serafol® foil	803030	100
SIFIN foil	BG 1713	100
Serafol® sticks	803020	50
SIFIN sticks	BG 1712	100

Test material

Capillary-, vein-, or stored donation blood (segment)

Test procedure

Method 1



1. Open foil wrapper and take out twin card.
2. When using as single card, fold twin card along perforation and tear apart prior to use.
3. Fill in recipient and blood donation data.
4. Add drops of material to be tested on reaction fields.

Recipient

Add one drop of recipient blood (50 µl) with a pipette or directly from the finger tip to each reaction field of one row (see fig. 1 and 2). Add one drop isotonic saline solution or water to each reaction field.

Donor (segment)

Add one drop of isotonic saline solution or water and then one drop (50 µl) donor blood to each reaction field in one row.



Auto control

The field "Blood" in each row is used to perform an auto control.

Add one drop of isotonic saline solution or water plus one drop recipient or donor blood to the field (see also instructions "recipient" or "donor").

5. Continue with test procedure immediately. Do not let blood dry.
6. Mix each field with the mixing stick for approx. 30 seconds until the reagent is completely dissolved. While mixing, spread material to be tested over the entire reaction field (see fig. 3). Use a new mixing stick for each reaction field or carefully clean the mixing stick in order to avoid antibody spreading. Use this procedure also for the auto control.

7. Keep the card in circular motion for approx. 30 to 60 seconds so that the drops circle in the reaction field (see fig. 4).

8. Observe for agglutination on the test fields and write test result on card.

Test results:

- Agglutination = positive reaction
- No agglutination = negative reaction

The eight possible blood type combinations (A, B, AB and O with positive or negative D rhesus characteristic) are listed in the following chart:

Anti-A	Anti-B	Anti-D	Auto control	Blood group
+	-	+	-	A Rh pos.
+	-	-	-	A Rh neg.
-	+	+	-	B Rh pos.
-	+	-	-	B Rh neg.
+	+	+	-	AB Rh pos.
+	+	-	-	AB Rh neg.
-	-	+	-	O Rh pos.
-	-	-	-	O Rh neg.
+	+	+	+	Test void

If the Serafol card result does not match the recipient or blood donation **DO NOT TRANSFUSE**, immediately find the cause and contact a transfusion medicine specialist!

To rely on the test result, a proper auto control is required. If the auto control field shows agglutination, the test result is not valid and must be re-checked.

Method 2:

1. to 3. see **method 1**
4. Add 1 drop of isotonic saline solution or water to the reaction fields.
5. Add 1 large drop (>80 µl) recipient or donor blood to one of the fields labelled "blood" on the card.
6. Dip a mixing stick as flat as possible into the blood drop and turn the stick until a ball of blood adheres to the end of the stick.
7. Dip the stick end with the blood into the drop of isotonic saline solution or water on the reaction field with Anti-A and mix until the reagent is completely dissolved (approx. 15 seconds). While mixing, spread material to be tested over the entire reaction field. Repeat this procedure with a new mixing stick on the reaction field Anti-B (or Anti-B and Anti-D) in the same row.

8. Set card aside for approx. 30 seconds.
9. Cautiously wave the card or turn to all sides. After 60 to 90 seconds observe for agglutination and write results on card.

Documentation

- Following transfusion the data should be filled in and the treating physician should sign.
- For filing and documentation, the reaction mixes can be dried and covered by a self-adhesive, see-through foil cover. During the drying the cards should be stored horizontally and should not be moved.

Notes

- Manual techniques are to be performed according to the manufacturer's instructions. Each deviation from these instructions is the sole responsibility of the user.
- Used test must be discarded as hazardous material.
- The bedside test is no proper substitution for complete blood typing or serological compatibility tests (cross matching).
- Because of the reagents' biotechnological production the risk of contamination with infectious agents is almost impossible. As some content derives from animal material (fetal calf serum, stabilizer) all test reagents must be regarded as potential transmitters of infectious diseases and be handled accordingly.
- As the test reagents contain sodium azide avoid any contact with the skin or the mucosa!

Limits of the procedure

- Due to the use of monoclonal reagents usually red cells with weak A-characteristics are detected: A₃ bloods with normal to weak reactions, A_x bloods weak positive to negative.
- The Anti-D reagent is a human monoclonal antibody of immunoglobulin-isotype M (clone BS226). It does not detect category VI cells.
- D^{weak}-bloods react weak or negative with the Anti-D reagent on Serafol. If reactions with the Anti-D reagent are negative or weak positive, this may indicate that the tested blood is D^{weak} or another rare variation. These antigens are detected during laboratory testing and have to be taken into account.

Sources of errors

Causes for false positive results

- Blood samples may occasionally react by rouleaux formation¹, which can be mistaken for a weak agglutination and may incorrectly be read as a positive result. This phenomenon has non-immunological causes and usually occurs when blood is directly and without isotonic saline solution or water added to the reaction fields and mixed. Rouleaux formation also occurs in heparin blood, blood from patients treated with plasma expanders (e.g. dextran or hydroxyethyl starch) as well as in blood from patients with plasmacytoma (high protein content, changes in protein content), oncological disease (abnormal blood count) or coagulation dysfunctions.

For testing these patients, first 1 drop of isotonic saline solution and then a drop of blood of the same size should be added to the reaction fields.

If rouleaux formation is suspected as the cause for false positive results or other reaction results for either the recipient or donor (auto control is not clear), the test should be repeated with method 2, which usually avoids this phenomenon.

¹ Reference: Issitt, Peter, D.: Applied Blood Group Serology, 4th Edition 1998, page 1134

- Drying out may be incorrectly interpreted as a positive reaction. *Perform test without delays.*
- Cold-reacting autoantibodies with high temperature amplitude may also cause a false positive result. These antibodies are detected during laboratory testing and have to be taken into account.

Causes for false negative results

- Blood drop too small or too large. *Repeat test precisely following instructions.*
- With low red cell concentrations (lower than 40 %) in the donor- or recipient blood drop, method 2 may result in delayed or weak reactions. *Repeat with method 1.*
- Expiry date has passed or card was stored improperly (over 30 °C for a longer time period). *Repeat test with new card.*
- Test was performed with red cell suspension. *Repeat test with vein-, capillary-, or donation blood.*
- Blood with a hematocrit under 15 % may cause weak or false negative reactions.

Package size

Serafol® ABO:	Single use kit (1 twin card, 2 cover foils, 2 mixing sticks, 1 reaction table)	Catalogue no.: 803100
		BG 1721
Serafol® ABO:	50 twin cards	Catalogue no.: 803110
		BG 1722
Serafol® ABO+D:	50 twin cards	Catalogue no.: 803120
		BG 1723

Symbol Explanation:

	Caution, refer to accompanying documents		Use by YYYY-MM (MM = end of month)
	Batch Code (Lot)		Upper limit of temperature
	Catalogue number		In Vitro Diagnostic Medical Device

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