

UREA-B



Urea.Berthelot. Enzymatic colorimetric
Quantitative Determination of urea
 Only for *in vitro* use in clinical laboratory
 Store at 2-8°C


Ref.:URB-030

2 x 125 mL

PRINCIPLE OF THE METHOD

Urea in the sample is hydrolyzed enzymatically into ammonia (NH₄⁺) and carbon dioxide (CO₂).

Ammonia ions formed reacts with salicylate and hypochlorite (NaClO), in presence of the catalyst nitroprusside, to form a green indophenol:



The intensity of the color formed is proportional to the urea concentration in the sample^{1,2,3}.

CLINICAL SIGNIFICANCE

Urea is the final result of the metabolism of proteins; it is formed in the liver from its destruction.

Elevated urea can appear in blood (uremia) in: diets with excess of proteins, renal diseases, heart failure, gastrointestinal hemorrhage, dehydration or renal obstruction^{1,6,7}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

| | | |
|-----------------|-------------------------------|------------|
| R 1 | Phosphate pH 6.7 | 50 mmol/L |
| Buffer | EDTA | 2 mmol/L |
| | Sodium salicylate | 400 mmol/L |
| | Sodium nitroprusside | 10 mmol/L |
| R 2 | Sodium hypochlorite (NaClO) | 140 mmol/L |
| NaClO | Sodium hydroxide | 150 mmol/L |
| R 3 | Urease | 30000 U/L |
| Enzymes | | |
| UREA CAL | Urea aqueous primary standard | 50 mg/dL |

PRECAUTIONS

R1 (sodium nitroprusside): **Xn, N: Harmful. Environmentally dangerous**

R26/27/28 Very toxic by inhalation, in contact with skin and if swallowed.

R32 Contact with acids liberates very toxic gas. R50 Harmful to aquatic organisms. R53 may cause long-term adverse effects in the aquatic environment. S7 Keep container tightly closed. S28 After contact with skin, wash immediately with plenty of water. S29 Do not empty into drains.

S60 This material and its container must be disposed of as hazardous waste.

S61 Avoid release to the environment. Refer to special instructions/safety data sheets.

R2 (Sodium hydroxide): **Xi, Irritant**. S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. S37/39 Wear suitable gloves and eye/face protection.

S45 In case of accident or if you feel unwell, sep medical advice immediately.

PREPARATION

- Working reagent (WR): Dissolve (→) one tablet R 3 Enzymes in one bottle of R 1 Buffer. Cap and mix gently to dissolve contents.

Stability: 4 weeks in the refrigerator (2-8°C) or 7 days at room temperature (15-25°C).

- R 2 NaClO is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

UREA CAL

Once open is stable up to 1 month when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm \geq 0.32.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 580 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment ^(Note 1).

SAMPLES

- Serum or heparinized plasma¹: Do not use ammonium salts or fluoride as anticoagulants.

- Urine¹: Dilute sample 1/50 in distilled water. Mix. Multiply results by 50 (dilution factor). Preserve urine samples at pH < 4.

Urea is stable at 2-8°C for 5 days;

PROCEDURE

- Assay conditions:
 Wavelength: 580 nm
 Cuvette: 1 cm light path
 Temperature: 37°C / 15-25°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:

| | Blank | Standard | Sample |
|-------------------------------------|-------|----------|--------|
| WR (mL) | 1.0 | 1.0 | 1.0 |
| Standard ^(Note 2-3) (□L) | -- | 10 | -- |
| Sample (□L) | -- | -- | 10 |
- Mix and incubate 5 min at 37°C or 10 min at room temperature (15-25°C).
- Pipette:

| | Blank | Standard | Sample |
|----------|-------|----------|--------|
| R 2 (mL) | 1.0 | 1.0 | 1.0 |
- Mix and incubate 5 min at 37°C or 10 min at room temperature (15-25°C).
- Read the absorbance (A) of the samples and calibrator, against the Blank. The colour is stable for at least 30 minutes at 15-25°C.

CALCULATIONS

$\frac{(A)_{\text{Sample}}}{(A)_{\text{Standard}}} \times 50$ (Standard conc.) = mg/dL urea in the sample

10 mg/L urea BUN divided by 0.466 = 21 mg/L urea = 0.36 mmol/L urea¹.

Conversion factor: mg/dL x 0.1665 = mmol/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures.

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

Serum : 15- 45 mg/dL (2.49-7.49 mmol/L)

Urine : 20 - 35 gr/24 h.

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0.3 mg/dL to linearity limit of 200 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

| | Intra-assay (n=20) | | Inter-assay (n=20) | |
|--------------|--------------------|------|--------------------|------|
| Mean (mg/dL) | 40.0 | 139 | 40.0 | 142 |
| SD | 1.27 | 3.50 | 1.86 | 3.75 |
| CV (%) | 3.17 | 2.50 | 4.64 | 2.63 |

Sensitivity: 1 mg/dL = 0.00505 A.

Accuracy: Results obtained using BSM reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0.9941.

Regression equation: y = 0.9972x + 0.011.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

It is recommended to use heparin as anticoagulant. Do not use ammonium salts or fluoride¹.

A list of drugs and other interfering substances with urea determination has been reported by Young et. al^{4,5}.

NOTES

- Glassware and distilled water must be free of ammonia and ammonium salts¹.
- Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
- Use clean disposable pipette tips for its dispensation.
- BSM has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

- Kaplan A. Urea. Kaplan A et al. *Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton* 1984; 1257-1260 and 437 and 418.
- Tabacco A et al. *Cin Chem* 1979; 25: 336-337.
- Fawcett J K et al. *J Clin Path* 1960; 13: 156-169.
- Young DS. *Effects of drugs on Clinical Lab. Tests*, 4th ed AACC Press, 1995.
- Young DS. *Effects of disease on Clinical Lab. Tests*, 4th ed AACC 2001.
- Burtis A et al. *Tietz Textbook of Clinical Chemistry*, 3rd ed AACC 1999.
- Tietz N W et al. *Clinical Guide to Laboratory Tests*, 3rd ed AACC 1995.